

Sublative Information

HOW DOES IT WORK?

The Sublative applicator delivers bipolar radio frequency energy via a matrix of 64 electrodes in the disposable applicator tip. The Radio Frequency energy generates fractional deep dermal heating in the area of the electrode matrix to prompt skin injury, thus causing a wound healing reaction. Mild to moderate wrinkles and superficial skin discoloration are reduced, and skin texture becomes more smooth and elastic.

IS IT A LASER?

It is neither laser- nor light-based. Instead, it harnesses the power of fractionated bi-polar radiofrequency (also known as gentle energy currents) to deeply penetrate the skin and kindle collagen production beneath the surface while keeping the outer layer of skin intact.

WHAT DOES IT TREAT?

Sublative results are immediate and continue to improve over the course of a few months:

- Improved elasticity
- Improved radiance
- More lifted
- More even tone
- Smoother, more supple texture
- Fewer acne scars, discolorations and sun damage

WHO SHOULD BE TREATED?

This treatment is a safe and effective solution for all skin types.

WHEN CAN PATIENTS EXPECT TO SEE RESULTS?

Post treatment, patients may notice tone and texture enhancements as a result of increased levels of beneficial collagen and new healthy skin cells which are produced. With the treatment's dermal impact, mild to moderate wrinkles and superficial skin discolorations are reduced. Typical treatment protocol is 3 treatments, 4 - 6 weeks apart. The total number of required treatment sessions depends on the skin condition of the patient.

HOW IS THE TREATMENT TOLERATED BY PATIENTS?

During treatment, many patients feel a warm, prickly sensation as energy enters their skin. After treatment a pink or red "sunburn" appearance and feeling is also common. Topical anesthetic ointments and cooling are often used to reduce discomfort during treatment.